



UNIVERSITY *of* MARYLAND
SCHOOL OF PHARMACY

Reclining languorously on my Phase Longue




It makes FT-ICR so
much more absorbing!

Phasing WTF?

- Phase corrected Absorption mode spectrum view offers an alternative to magnitude mode
- It's a post processing trick you can play on FT-ICR MS data which brings some major benefits
- Benefits
- How and Why
- New Advances

Benefits of absorption mode

- Improved resolving power
- Sample
 - Pyrolysis fluid from waste organic matter feedstock




Benefits of absorption mode

- Improved resolving power
- Sample
 - Pyrolysis fluid from waste organic matter feedstock




Been huntin' those peaks since 1850, dammit!



Benefits of absorption mode


- Improved resolving power
- Same improvements also in protein spectra
 - More peaks detectable



Protein coverage

- Improved signal to noise ratio
- Means
 - More fragments detected
 - More confident assignment
 - Better sequence coverage
 - More PTM information

Magnitude mode




- Cannot determine the charge state for these peak clusters

1439 1440 1441 1442 1443 1444 1445 1446 m/z

Protein coverage

- Improved signal to noise ratio
- Means
 - More fragments detected
 - More confident assignment
 - Better sequence coverage
 - More PTM information

Magnitude mode



Absorption mode

Improved Mass Accuracy

- Sample




- Talisker Malt Whisky
- 12T Solarix, + mode ESI

Peak Assignment by KMD algorithm^{1,2}

- Assigned peak mass accuracy

- Magnitude mode
0.080ppm \pm 0.104
- Absorption mode
0.066ppm \pm 0.084



1. Kilgour, Mackay, Langridge-Smith, O'Connor. Anal Chem, 84(17), 7431-7435

2. Kilgour, Mackay, Langridge-Smith, O'Connor. Anal Chem, 84(17), 7436-7439


Benefits of absorption mode

- Features
 - Improved resolving power
 - Improved signal-to-noise
 - Better mass accuracy
- Benefits
 - You will be able to see more peaks
 - You can assign with more confidence

– **It's free!!**





Phase shifts FT-ICR MS



Phase shifts FT-ICR MS



- Excite and delay cause phase shift



- By the beginning of detection
 - Ions with different frequencies have different starting angles
- This causes the peak distortion

Phase shifts FT-ICR MS



- Excite and delay cause phase shift



- By the beginning of detection
 - Ions with different frequencies have different starting angles
- This causes the peak distortion

Phase shifts FT-ICR MS


- Excite and delay cause phase shift




- By the beginning of detection
 - Ions with different frequencies have different starting angles
- This causes the peak distortion

Phase shifts FT-ICR MS

- Excite and delay cause phase shift




- Phase correction aim
 - Remove all accumulated phase




Phase shifts FT-ICR MS

- Excite and delay cause phase shift




- Phase correction aim
 - Remove all accumulated phase



Makes it look like all ions started at the same place

The fly in the ointment

- Magnitude mode spectrum
- Phase corrected absorption mode
- Un-phase corrected absorption mode



How do you make an absorption mode spectrum

- Record ion transient
- FFT
 - Extract magnitude and phase angle of all ion frequencies
- Calculate phase angle correction function
 - Quadratic relationship
 - $\varphi_f = Af^2 + Bf + C$

f – frequency
φ_f – phase angle correction
A, B & C – coefficients that are optimised
- Plot phase corrected absorption mode spectrum


Autophaser Method

- Simple method of producing Absorption mode spectra

- PCF

- Phase Correction Function

- $\varphi_f = Af^2 + Bf + C$




Autophaser Method

- Simple method of producing Absorption mode spectra

- Basis Phase Correction Function


- Estimate from spectrum*
- Load for previous spectrum*
- Generate from excite function

- *Also options in Bruker software




Autophaser Method

- Simple method of producing Absorption mode spectra



- Optimise Phase Correction Function for current spectrum
 - Account for ion population related frequency shifts




So, What's new?

- The perfidious baseline correction
 - And how to avoid it
- Non-quadratic phase correction functions
 - When are they needed
 - What causes them


Baseline Correction

Raw
absorption
mode




Baseline Correction

Raw
absorption
mode




Baseline Correction

Raw
absorption
mode




Baseline
corrected




Baseline Correction

Raw
absorption
mode

Magnitude
mode




Baseline
corrected



Relative peak height


- No baseline correction
- Detected baseline
- After baseline correction



Very difficult to produce accurate baseline fit fully automatically, reliably

Relative peak height

- No baseline correction
- Detected baseline
- After baseline correction




3s per pixel on a 50k pixel image = 42 extra hours of processing!



Takes ~3s (1M transient) – ~20 secs (8M transient) to compute


Relative peak height

- Magnitude mode
 - Hanning
 - As ref
- Absorption mode peaks
 - Before Corr
 - After Corr



Relative peak height


- Magnitude mode
 - Hanning
 - As ref
- Absorption mode peaks
 - Before Corr
 - After Corr



Baseline correction may perturb isotopic relative abundances


Relative peak height

- No baseline correction
 - Mean 0.98
 - SD 0.03



Relative peak height



- No baseline correction
 - Mean 0.98
 - SD 0.03
- After baseline correction
 - Mean 1.02
 - SD 0.05



Baseline correction may cause more problems than it solves



Apodization games

- Fourier transform
 - Assumes the waveform is infinitely long
 - Generates spectral leakage on finite signals
- Apodization or windowing
 - Helps limit spectral leakage
- Half-Hanning function
 - Commonly used for Absorption mode spectra





Apodization games

- Fourier transform
 - Assumes the waveform is infinitely long
 - Generates spectral leakage on finite signals
- Apodization or windowing
 - Helps limit spectral leakage
- Hanning function
 - Commonly used for Magnitude mode spectra




Results

- Half-Hanning apodized absorption mode
- Full Hanning apodized absorption mode
- Baseline is now flat
 - On the wide scale





Results

- Half-Hanning apodized absorption mode
- Full Hanning apodized absorption mode
- Baseline is now flat
 - On the wide scale
 - Also on zoom






Results

- Half-Hanning apodized absorption mode
- Full Hanning apodized absorption mode
- Baseline is now flat
 - On the wide scale
 - Also on zoom



Results

- Half-Hanning apodized absorption mode
- Full Hanning apodized absorption mode
- Concern






C'mon bruv!

MS don't have
no negative
side lobes!

Results




- Half-Hanning apodized absorption mode
- Full Hanning apodized absorption mode
- Concern



No no,
homeboy. The
sidelobes are
true, innit?

Results


- Half-Hanning apodized absorption mode
- Full Hanning apodized absorption mode
- Concern




But you can use an adaptive filter to reduce them

Results

- Half-Hanning apodized absorption mode
- Full Hanning apodized absorption mode
- Magnitude Mode





Gosh – that's much better than magnitude mode!




Results

- Relative peak intensity
 - Absorption mode vs Magnitude mode
 - Mean – 1.00
 - SD – 0.01





Benefits

- Avoiding baseline correction
 - Speeds up processing
 - Improves relative abundance accuracy
- Should be very beneficial in FT-ICR imaging
 - Also makes peak detection easier



Phase shifts FT-ICR MS



- Excite and delay cause phase shift



- By the beginning of detection
 - Ions with different frequencies have different starting angles
- This causes the peak distortion

Phase shifts FT-ICR MS



- Excite and delay cause phase shift



- By the beginning of detection
 - Ions with different frequencies have different starting angles
- This causes the peak distortion

Phase shifts FT-ICR MS


- Excite and delay cause phase shift




- By the beginning of detection
 - Ions with different frequencies have different starting angles
- This causes the peak distortion

Phase shifts FT-ICR MS

- Excite and delay cause phase shift




- Phase correction aim
 - Remove all accumulated phase




Phase shifts FT-ICR MS

- Excite and delay cause phase shift



- Phase correction aim
 - Remove all accumulated phase




Makes it look like all ions started at the same place

Source of phase shifts

- Frequency sweep and delay times

$$\varphi = \frac{-\omega^2}{R} + \omega \left(\frac{\omega_{Final}}{R} + t_{Delay} \right)$$



- Legend

- ϕ_i – initial phase angle (rad) of excite
- ω_x – detected ion frequency (rad.s^{-1})
- ω_0 – initial frequency (rad.s^{-1}) of excite
- ω_{Final} – final frequency (rad.s^{-1}) of excite
- R – excitation sweep rate
- t_{Delay} – delay time (s)


Source of phase shifts

- Initial phase

$$\varphi_t = \frac{\omega_t^2}{2R} + \frac{-\omega_0^2}{R} + \varphi_i$$

- Frequency sweep and delay times

$$\varphi = \frac{-\omega^2}{R} + \omega \left(\frac{\omega_{Final}}{R} + t_{Delay} \right)$$



- Legend

- ϕ_i – initial phase angle (rad) of excite
- ω_x – detected ion frequency (rad.s⁻¹)
- ω_0 – initial frequency (rad.s⁻¹) of excite
- ω_{Final} – final frequency (rad.s⁻¹) of excite
- R – excitation sweep rate
- t_{Delay} – delay time (s)

Source of phase shifts

- Initial phase

$$\varphi_t = \frac{\omega_t^2}{2R} + \frac{-\omega_0^2}{R} + \varphi_i$$

- Frequency sweep and delay times

$$\varphi = \frac{-\omega^2}{R} + \omega \left(\frac{\omega_{Final}}{R} + t_{Delay} \right)$$

- Legend

- ϕ_i – initial phase angle (rad) of excite
- ω_x – detected ion frequency (rad.s⁻¹)
- ω_0 – initial frequency (rad.s⁻¹) of excite
- ω_{Final} – final frequency (rad.s⁻¹) of excite
- R – excitation sweep rate
- t_{Delay} – delay time (s)

$$\varphi_x = \frac{-\omega_x^2}{2R} + \omega_x \left(\frac{\omega_{Final}}{R} + t_{Delay} \right) + \frac{-\omega_0^2}{R} + \varphi_i$$

Source of phase shifts

- Initial phase

$$\varphi_t = \frac{\omega_t^2}{2R} + \frac{-\omega_0^2}{R} + \varphi_i$$

- Frequency sweep and delay times

$$\varphi = \frac{-\omega^2}{R} + \omega \left(\frac{\omega_{Final}}{R} + t_{Delay} \right)$$

- Total equation

$$\varphi_x = \frac{-\omega_x^2}{2R} + \omega_x \left(\frac{\omega_{Final} - \Delta\omega_K}{R} + t_{Delay} \right) + \left(\frac{-\omega_0^2}{2R} + \varphi_i \right) + \frac{\Delta\omega_K \omega_{Final}}{R}$$

You never really know this, accurately

- Complication

- Frequency shifts

- $\bar{\omega}_i = \omega_x + \Delta\omega_K$

- Legend

- φ_i – initial phase angle (rad) of excite
- ω_x – detected ion frequency (rad.s⁻¹)
- ω_0 – initial frequency (rad.s⁻¹) of excite
- ω_{Final} – final frequency (rad.s⁻¹) of excite
- R – excitation sweep rate
- t_{Delay} – delay time (s)

Source of phase shifts

- Initial phase

$$\varphi_t = \frac{\omega_t^2}{2R} + \frac{-\omega_0^2}{R} + \varphi_i$$

- Frequency sweep and delay times

$$\varphi = \frac{-\omega^2}{R} + \omega \left(\frac{\omega_{Final}}{R} + t_{Delay} \right)$$

- Total equation

$$\varphi_x = \frac{-\omega_x^2}{2R} + \omega_x \left(\frac{\omega_{Final} - \Delta\omega_K}{R} + t_{Delay} \right) + \left(\frac{-\omega_0^2}{2R} + \varphi_i \right) + \frac{\Delta\omega_K \omega_{Final}}{R}$$

You never really know this, accurately

- Complication

- Frequency shifts

- $\bar{\omega}_i = \omega_x + \Delta\omega_K$

- Legend

- φ_i – initial phase angle (rad) of excite
- ω_x – detected ion frequency (rad.s⁻¹)
- ω_0 – initial frequency (rad.s⁻¹) of excite
- ω_{Final} – final frequency (rad.s⁻¹) of excite
- R – excitation sweep rate
- t_{Delay} – delay time (s)


Phase correction by fitting

- Full quadratic fitting is required
 - To allow for $\Delta\omega_k$
- Phase correction function
 - $\varphi_f = Af^2 + Bf + C$
- Can solve for A , B and C (if you know pairs of φ and f)
- With no phase correction
 - $Re = Mag * \cos \theta$
- With phase correction
 - $Re_f = Mag_f \times \cos(\theta_f - \varphi_f)$
 - φ_f is the phase correction factor

An unexpected pleasure

- The file that wouldn't phase correct ...
- Residual phase error is not quadratic
 - Original phase error was of a higher order
 - Probably third order
- (In this case caused by a non-linear frequency sweep)


Residual phase error after optimization



An unexpected pleasure



- The file that wouldn't phase correct ...
- Residual phase error is not quadratic
 - Original phase error was of a higher order
 - Probably third order
- (In this case caused by a non-linear frequency sweep)

Residual phase error after optimization



Solution

- Adapt Autophaser to solve n^{th} order PCFs
 - Simple adaptation to GA to fit curve to $(n+1)$ nodes instead of 3



- 3rd order fitting actually makes it easier to phase correct other files too

Conclusions

- An new apodization approach
 - Removes the necessity of baseline correction
 - Better abundance sensitivity
 - Improved processing times
 - Will be key for absorption imagery
- n^{th} order PCFs
 - Can be solved!
 - May be caused by unusual excitation functions



is more fun even
than riding my
Vespa down a
cobble street!

